

**Clinical Agreement Study for
SARS-CoV-2 Neutralizing Antibody
Fast Test Kit (Immunofluorescence Assay)**

CONTENTS

| | |
|---|----|
| 1. Purpose..... | 3 |
| 2. Experimental Materials | 3 |
| 2.1 Trial reagent..... | 3 |
| 2.2 Comparator reagent | 3 |
| 3. Experimental sample | 3 |
| 3.1 Sample source and amount..... | 3 |
| 3.2 Sample information | 4 |
| 4. Test process | 16 |
| 5. Statistical methods | 16 |
| 5.1 Calculation method of positive/negative percent agreement | 16 |
| 5.2 Calculation method of 95% confidence interval | 17 |
| 6. Test results..... | 17 |
| 6.1 Clinical agreement study result | 17 |
| 7. Statistical analysis..... | 26 |
| 7.1 Clinical agreement study results analysis..... | 26 |
| 8. Conclusion | 27 |

1. Purpose

SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) is developed by Getein Biotech, Inc. The test is intended for the qualitative detection of anti-SARS-CoV-2 neutralizing antibodies in serum, plasma or whole blood samples. SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection and should not be used to diagnose acute SARS-CoV-2 infection.

The aim of the clinical agreement study was to compare and evaluate the clinical performance of SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) with EUA-authorized tests.

2. Experimental Materials

2.1 Trial reagent

Name: SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay)

Specification: 25 tests per box

Lot no.: YXGN210001W (Manufacturing date: January 06, 2021)

Manufacturer: Getein Biotech, Inc.

2.2 Comparator reagent

Name: cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit (EUA Approved)

Specifications: 96 Tests

Manufacturer: GenScript Biotech Co., Ltd.

Lot no.: A201209 (Expiration date: December 20, 2021)

3. Experimental sample

3.1 Sample source and amount

A total of 409 samples were tested using SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay). The combined cohort was consisted of healthy people (n=257, RT-PCR negative) and volunteers who have been immunized by injecting two doses of vaccine against SARS-CoV-2 (n=152). Serum samples of vaccine group were collected at 28 days after the second dose of vaccine.

Samples of serum were collected. The healthy people samples were consented and collected from Nanjing Jiangbei People's Hospital and the volunteer's samples by injecting vaccine were collected from Getein Biotech, Inc. and Nanjing Jiangbei People's Hospital based on the sample inclusion criteria. The studies were performed at January 2021.

3.2 Sample information

The samples were tested with GenScript Biotech Co., Ltd's cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit for SARS-CoV-2 neutralization antibody. In this study, the GenScript Biotech Co., Ltd testing positive samples were considered as positive, the testing negative samples were considered as negative.

3.2.1 Samples of Clinical agreement study

A total of 409 samples were used for the Clinical Agreement Study. The detailed sample information is as the following:

The detailed information of 257 healthy people samples data and sources are listed below in the Table 1.

Table 1 The 257 healthy people samples information of clinical agreement study from the Nanjing Jiangbei People's Hospital

| No. | Sample ID | Gender | Age | Sample type | Sample source |
|-----|-----------|--------|-----|-------------|------------------------------------|
| 1 | S470025 | Male | 17 | Serum | Nanjing Jiangbei People's Hospital |
| 2 | S470008 | Female | 14 | | |
| 3 | S473257 | Male | 57 | | |
| 4 | S476254 | Male | 79 | | |
| 5 | S470125 | Female | 24 | | |
| 6 | S476354 | Male | 57 | | |
| 7 | S471025 | Female | 75 | | |
| 8 | S473156 | Female | 48 | | |
| 9 | S472104 | Male | 66 | | |
| 10 | S472531 | Female | 78 | | |
| 11 | S470075 | Male | 14 | | |
| 12 | S471034 | Female | 20 | | |
| 13 | S471124 | Male | 41 | | |
| 14 | S470313 | Female | 8 | | |
| 15 | S477739 | Male | 33 | | |
| 16 | S471641 | Female | 76 | | |
| 17 | S476291 | Male | 34 | | |
| 18 | S477613 | Male | 7 | | |
| 19 | S477547 | Male | 9 | | |
| 20 | S472570 | Female | 7 | | |
| 21 | S477639 | Male | 37 | | |
| 22 | S478189 | Male | 43 | | |

| | | | |
|----|---------|--------|----|
| 23 | S475576 | Female | 11 |
| 24 | S476956 | Male | 32 |
| 25 | S477568 | Male | 12 |
| 26 | S476940 | Female | 58 |
| 27 | S471619 | Female | 92 |
| 28 | S476493 | Male | 43 |
| 29 | S471454 | Male | 52 |
| 30 | S474205 | Male | 26 |
| 31 | S470007 | Female | 41 |
| 32 | S479521 | Male | 24 |
| 33 | S472826 | Female | 34 |
| 34 | S473021 | Male | 8 |
| 35 | S474966 | Male | 12 |
| 36 | S475023 | Male | 56 |
| 37 | S475722 | Female | 54 |
| 38 | S472016 | Male | 84 |
| 39 | S477837 | Female | 75 |
| 40 | S472647 | Female | 14 |
| 41 | S472163 | Male | 58 |
| 42 | S474983 | Male | 42 |
| 43 | S479566 | Female | 79 |
| 44 | S479141 | Male | 39 |
| 45 | S476713 | Female | 21 |
| 46 | S473363 | Male | 21 |
| 47 | S473338 | Female | 77 |
| 48 | S471304 | Male | 24 |
| 49 | S477205 | Male | 15 |
| 50 | S478990 | Female | 82 |
| 51 | S473027 | Male | 76 |
| 52 | S479358 | Female | 19 |
| 53 | S476855 | Female | 29 |
| 54 | S474914 | Male | 34 |
| 55 | S477793 | Male | 31 |
| 56 | S476871 | Male | 52 |

| | | | |
|----|---------|--------|----|
| 57 | S477591 | Female | 30 |
| 58 | S477625 | Male | 75 |
| 59 | S478298 | Female | 14 |
| 60 | S476990 | Male | 68 |
| 61 | S473915 | Male | 53 |
| 62 | S475388 | Female | 72 |
| 63 | S479392 | Male | 79 |
| 64 | S472070 | Male | 14 |
| 65 | S479543 | Male | 53 |
| 66 | S476304 | Female | 23 |
| 67 | S478389 | Male | 77 |
| 68 | S473500 | Female | 38 |
| 69 | S475104 | Male | 41 |
| 70 | S475595 | Female | 48 |
| 71 | S475617 | Male | 74 |
| 72 | S476163 | Male | 34 |
| 73 | S474991 | Female | 51 |
| 74 | S478618 | Male | 16 |
| 75 | S478292 | Male | 57 |
| 76 | S473196 | Female | 64 |
| 77 | S477450 | Male | 73 |
| 78 | S477701 | Male | 14 |
| 79 | S479113 | Male | 61 |
| 80 | S472751 | Female | 43 |
| 81 | S477686 | Female | 49 |
| 82 | S479780 | Male | 69 |
| 83 | S472855 | Male | 27 |
| 84 | S470629 | Male | 20 |
| 85 | S479106 | Female | 51 |
| 86 | S479186 | Male | 36 |
| 87 | S475911 | Male | 32 |
| 88 | S479254 | Female | 62 |
| 89 | S471417 | Female | 64 |
| 90 | S470609 | Male | 30 |

| | | | |
|-----|---------|--------|----|
| 91 | S470273 | Male | 36 |
| 92 | S470785 | Male | 27 |
| 93 | S472317 | Female | 49 |
| 94 | S479124 | Male | 33 |
| 95 | S474178 | Male | 75 |
| 96 | S478444 | Female | 63 |
| 97 | S471263 | Female | 69 |
| 98 | S470143 | Male | 33 |
| 99 | S470447 | Male | 27 |
| 100 | S471040 | Male | 62 |
| 101 | S474369 | Female | 80 |
| 102 | S475657 | Female | 15 |
| 103 | S477997 | Male | 54 |
| 104 | S470290 | Male | 35 |
| 105 | S476907 | Female | 76 |
| 106 | S473240 | Male | 55 |
| 107 | S474737 | Male | 36 |
| 108 | S473218 | Female | 82 |
| 109 | S478393 | Male | 65 |
| 110 | S479487 | Female | 77 |
| 111 | S472566 | Male | 81 |
| 112 | S473558 | Male | 48 |
| 113 | S478069 | Female | 37 |
| 114 | S477800 | Male | 22 |
| 115 | S476683 | Male | 18 |
| 116 | S475643 | Female | 19 |
| 117 | S472399 | Male | 42 |
| 118 | S470593 | Female | 68 |
| 119 | S474765 | Female | 39 |
| 120 | S477504 | Male | 20 |
| 121 | S473573 | Male | 36 |
| 122 | S473024 | Female | 59 |
| 123 | S479962 | Female | 46 |
| 124 | S479693 | Female | 75 |

| | | | |
|-----|---------|--------|----|
| 125 | S471024 | Male | 49 |
| 126 | S476201 | Male | 64 |
| 127 | S477744 | Male | 54 |
| 128 | S470056 | Male | 52 |
| 129 | S475161 | Male | 45 |
| 130 | S472397 | Male | 19 |
| 131 | S470870 | Female | 63 |
| 132 | S475873 | Female | 39 |
| 133 | S474936 | Female | 67 |
| 134 | S476410 | Male | 72 |
| 135 | S478448 | Male | 25 |
| 136 | S476554 | Male | 22 |
| 137 | S477258 | Female | 65 |
| 138 | S475946 | Male | 64 |
| 139 | S478071 | Male | 73 |
| 140 | S473364 | Female | 63 |
| 141 | S476119 | Male | 79 |
| 142 | S479002 | Female | 65 |
| 143 | S470807 | Male | 37 |
| 144 | S477336 | Male | 19 |
| 145 | S475990 | Female | 28 |
| 146 | S475687 | Female | 29 |
| 147 | S473752 | Male | 65 |
| 148 | S475451 | Male | 14 |
| 149 | S476955 | Male | 78 |
| 150 | S472230 | Female | 61 |
| 151 | S477513 | Female | 38 |
| 152 | S470387 | Male | 6 |
| 153 | S475938 | Male | 71 |
| 154 | S475500 | Female | 81 |
| 155 | S474886 | Female | 42 |
| 156 | S475042 | Male | 70 |
| 157 | S477940 | Male | 57 |
| 158 | S470355 | Female | 57 |

| | | | |
|-----|---------|--------|----|
| 159 | S470900 | Female | 79 |
| 160 | S477572 | Male | 62 |
| 161 | S473643 | Female | 68 |
| 162 | S478650 | Male | 19 |
| 163 | S476116 | Female | 26 |
| 164 | S475810 | Male | 67 |
| 165 | S476986 | Male | 71 |
| 166 | S477199 | Male | 79 |
| 167 | S477935 | Female | 18 |
| 168 | S479357 | Male | 47 |
| 169 | S473104 | Male | 59 |
| 170 | S472752 | Female | 27 |
| 171 | S475556 | Male | 63 |
| 172 | S478734 | Male | 15 |
| 173 | S472941 | Male | 66 |
| 174 | S472240 | Male | 60 |
| 175 | S475000 | Female | 27 |
| 176 | S477376 | Female | 68 |
| 177 | S475735 | Male | 62 |
| 178 | S473474 | Male | 22 |
| 179 | S476234 | Male | 17 |
| 180 | S470368 | Male | 54 |
| 181 | S477938 | Female | 6 |
| 182 | S472433 | Male | 42 |
| 183 | S477965 | Female | 65 |
| 184 | S479301 | Male | 54 |
| 185 | S476584 | Female | 64 |
| 186 | S475542 | Male | 72 |
| 187 | S479578 | Male | 11 |
| 188 | S479223 | Female | 69 |
| 189 | S479668 | Female | 74 |
| 190 | S479394 | Female | 43 |
| 191 | S478014 | Male | 20 |
| 192 | S475919 | Male | 15 |

| | | | |
|-----|---------|--------|----|
| 193 | S472559 | Male | 38 |
| 194 | S478549 | Female | 50 |
| 195 | S474195 | Male | 43 |
| 196 | S474360 | Female | 76 |
| 197 | S478008 | Male | 35 |
| 198 | S472337 | Male | 54 |
| 199 | S472157 | Female | 58 |
| 200 | S477991 | Male | 53 |
| 201 | S470583 | Female | 27 |
| 202 | S471125 | Male | 12 |
| 203 | S472829 | Female | 78 |
| 204 | S470403 | Male | 35 |
| 205 | S471872 | Female | 21 |
| 206 | S471044 | Male | 37 |
| 207 | S477221 | Female | 44 |
| 208 | S470621 | Female | 21 |
| 209 | S473087 | Male | 32 |
| 210 | S478046 | Female | 6 |
| 211 | S470750 | Male | 9 |
| 212 | S479133 | Female | 28 |
| 213 | S471601 | Male | 64 |
| 214 | S472510 | Female | 74 |
| 215 | S470598 | Female | 65 |
| 216 | S473639 | Male | 47 |
| 217 | S472975 | Male | 36 |
| 218 | S471681 | Male | 61 |
| 219 | S474568 | Female | 32 |
| 220 | S474417 | Female | 25 |
| 221 | S476386 | Male | 63 |
| 222 | S472151 | Male | 41 |
| 223 | S477558 | Female | 22 |
| 224 | S477759 | Male | 36 |
| 225 | S474535 | Female | 79 |
| 226 | S478397 | Female | 25 |

| | | | |
|-----|---------|--------|----|
| 227 | S473088 | Male | 12 |
| 228 | S470739 | Male | 36 |
| 229 | S470556 | Male | 41 |
| 230 | S475313 | Female | 68 |
| 231 | S476196 | Male | 60 |
| 232 | S471383 | Male | 56 |
| 233 | S471969 | Female | 75 |
| 234 | S476392 | Female | 48 |
| 235 | S478093 | Male | 15 |
| 236 | S475441 | Male | 10 |
| 237 | S477819 | Female | 15 |
| 238 | S476663 | Male | 37 |
| 239 | S474724 | Female | 44 |
| 240 | S470163 | Male | 24 |
| 241 | S473374 | Female | 76 |
| 242 | S472014 | Male | 20 |
| 243 | S473722 | Female | 60 |
| 244 | S470924 | Female | 72 |
| 245 | S478073 | Female | 60 |
| 246 | S476318 | Male | 12 |
| 247 | S479875 | Female | 33 |
| 248 | S477760 | Female | 47 |
| 249 | S471033 | Male | 71 |
| 250 | S479257 | Female | 18 |
| 251 | S473112 | Female | 60 |
| 252 | S471613 | Male | 38 |
| 253 | S473817 | Female | 31 |
| 254 | S475208 | Female | 59 |
| 255 | S470698 | Male | 62 |
| 256 | S473414 | Male | 10 |
| 257 | S471552 | Female | 21 |

The detailed information of 152 volunteers (who have been immunized by injecting two doses of vaccine against SARS-CoV-2) samples data and sources are listed below in the Table 2.

Table 2 The sample information of 152 volunteers' samples

| No. | Sample ID | Gender | Age | Sample type | Volunteer source |
|-----|-----------|--------|-----|-------------|------------------------------------|
| 1 | S513164 | Female | 36 | Serum | Nanjing Jiangbei People's Hospital |
| 2 | S550045 | Female | 47 | | |
| 3 | S511159 | Male | 49 | | |
| 4 | S514157 | Female | 20 | | |
| 5 | S535624 | Male | 44 | | |
| 6 | S542825 | Female | 49 | | |
| 7 | S518563 | Male | 22 | | |
| 8 | S539775 | Male | 35 | | |
| 9 | S524989 | Female | 29 | | |
| 10 | S503466 | Female | 23 | | |
| 11 | S553136 | Male | 39 | | |
| 12 | S549251 | Male | 27 | | |
| 13 | S547035 | Female | 31 | | |
| 14 | S536878 | Female | 45 | | |
| 15 | S547249 | Female | 44 | | |
| 16 | S506840 | Male | 28 | | |
| 17 | S507842 | Female | 36 | | |
| 18 | S566038 | Female | 30 | | |
| 19 | S506292 | Male | 33 | | |
| 20 | S557506 | Male | 46 | | |
| 21 | S547257 | Male | 46 | | |
| 22 | S501290 | Female | 21 | | |
| 23 | S534916 | Female | 36 | | |
| 24 | S555978 | Female | 30 | | |
| 25 | S531674 | Male | 34 | | |
| 26 | S521406 | Female | 45 | | |
| 27 | S522183 | Male | 37 | | |
| 28 | S557341 | Female | 49 | | |
| 29 | S557566 | Female | 50 | | |
| 30 | S502151 | Male | 46 | | |
| 31 | S564978 | Female | 33 | | |
| 32 | S527532 | Male | 37 | | |

| | | | |
|----|---------|--------|----|
| 33 | S567357 | Female | 40 |
| 34 | S569672 | Male | 45 |
| 35 | S557509 | Female | 22 |
| 36 | S512594 | Male | 49 |
| 37 | S552563 | Female | 28 |
| 38 | S559789 | Male | 20 |
| 39 | S530763 | Female | 36 |
| 40 | S511158 | Male | 24 |
| 41 | S539063 | Female | 41 |
| 42 | S524045 | Female | 21 |
| 43 | S563778 | Male | 23 |
| 44 | S558088 | Female | 21 |
| 45 | S514071 | Male | 38 |
| 46 | S562071 | Female | 48 |
| 47 | S533450 | Female | 26 |
| 48 | S560602 | Male | 21 |
| 49 | S539696 | Female | 34 |
| 50 | S513867 | Male | 29 |
| 51 | S506056 | Male | 44 |
| 52 | S567823 | Female | 34 |
| 53 | S532216 | Male | 25 |
| 54 | S541848 | Female | 22 |
| 55 | S565421 | Female | 45 |
| 56 | S531101 | Male | 20 |
| 57 | S509082 | Female | 28 |
| 58 | S520230 | Female | 33 |
| 59 | S505999 | Male | 48 |
| 60 | S542741 | Female | 38 |
| 61 | S564448 | Male | 30 |
| 62 | S500522 | Female | 49 |
| 63 | S554527 | Male | 41 |
| 64 | S551319 | Female | 20 |
| 65 | S508357 | Male | 25 |
| 66 | S560772 | Female | 21 |

| | | | |
|-----|---------|--------|----|
| 67 | S548614 | Male | 29 |
| 68 | S504628 | Female | 24 |
| 69 | S501627 | Male | 47 |
| 70 | S540914 | Female | 33 |
| 71 | S565909 | Male | 44 |
| 72 | S535245 | Male | 36 |
| 73 | S554291 | Male | 42 |
| 74 | S567922 | Female | 22 |
| 75 | S514891 | Female | 38 |
| 76 | S532124 | Male | 40 |
| 77 | S567144 | Male | 34 |
| 78 | S514023 | Male | 38 |
| 79 | S502824 | Female | 34 |
| 80 | S537643 | Male | 27 |
| 81 | S533744 | Male | 38 |
| 82 | S552009 | Female | 50 |
| 83 | S500910 | Male | 43 |
| 84 | S550233 | Female | 38 |
| 85 | S529420 | Female | 34 |
| 86 | S568011 | Male | 20 |
| 87 | S568154 | Male | 42 |
| 88 | S557460 | Male | 30 |
| 89 | S518685 | Female | 27 |
| 90 | S563651 | Female | 46 |
| 91 | S557285 | Female | 42 |
| 92 | S555810 | Female | 47 |
| 93 | S538188 | Female | 34 |
| 94 | S519000 | Male | 39 |
| 95 | S516304 | Male | 32 |
| 96 | S545710 | Female | 45 |
| 97 | S521281 | Female | 25 |
| 98 | S530823 | Female | 32 |
| 99 | S565364 | Male | 47 |
| 100 | S562458 | Male | 50 |

| | | | |
|-----|---------|--------|----|
| 101 | S555117 | Female | 48 |
| 102 | S542380 | Male | 44 |
| 103 | S511206 | Male | 48 |
| 104 | S557798 | Female | 49 |
| 105 | S552736 | Male | 41 |
| 106 | S533004 | Male | 25 |
| 107 | S511162 | Male | 43 |
| 108 | S510089 | Male | 42 |
| 109 | S532254 | Female | 49 |
| 110 | S506297 | Female | 50 |
| 111 | S557476 | Female | 34 |
| 112 | S551231 | Male | 31 |
| 113 | S520068 | Female | 36 |
| 114 | S552794 | Male | 38 |
| 115 | S545502 | Male | 27 |
| 116 | S557538 | Female | 34 |
| 117 | S507353 | Male | 35 |
| 118 | S530192 | Female | 43 |
| 119 | S523997 | Male | 43 |
| 120 | S530458 | Female | 37 |
| 121 | S523242 | Female | 26 |
| 122 | S550714 | Male | 33 |
| 123 | S527507 | Male | 45 |
| 124 | S510608 | Male | 23 |
| 125 | S555694 | Female | 34 |
| 126 | S564136 | Female | 22 |
| 127 | S555403 | Male | 37 |
| 128 | S539594 | Female | 47 |
| 129 | S534609 | Female | 32 |
| 130 | S500909 | Male | 36 |
| 131 | S556215 | Female | 44 |
| 132 | S520221 | Female | 24 |
| 133 | S520093 | Female | 48 |
| 134 | S553199 | Female | 44 |

| |
|--|
| |
|--|

Getein Biotech, Inc.

| | | | |
|-----|---------|--------|----|
| 135 | S526734 | Male | 26 |
| 136 | S501436 | Female | 27 |
| 137 | S511877 | Male | 50 |
| 138 | S516012 | Male | 29 |
| 139 | S540804 | Female | 27 |
| 140 | S511845 | Male | 20 |
| 141 | S510532 | Female | 22 |
| 142 | S527460 | Female | 39 |
| 143 | S535409 | Male | 34 |
| 144 | S541180 | Female | 38 |
| 145 | S505230 | Male | 43 |
| 146 | S567996 | Female | 44 |
| 147 | S509345 | Male | 47 |
| 148 | S500710 | Male | 21 |
| 149 | S516760 | Female | 27 |
| 150 | S543837 | Female | 34 |
| 151 | S577513 | Female | 38 |
| 152 | S564542 | Male | 47 |

4. Test process

The clinical agreement study of SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) was evaluated by testing a total of 409 clinical samples from individual patients: The study was a double-blind and randomization experiment. A is responsible for the collection and numbering of all samples, and performing the comparator reagent testing. B is responsible for performing the trial reagent only (Note: A and B should not reveal the testing results to each other so to ensure the study was double-blinded). C is responsible for reveal and compare the testing results from A and B. D is responsible for the statistical analysis of the data.

5. Statistical methods

5.1 Calculation method of positive, negative and total percent agreement

Table 3 The test results of trial reagent and comparator reagent

| | |
|------------|--------------------|
| Experiment | comparator reagent |
|------------|--------------------|

| | | | |
|---------------|----------|----------|----------|
| | | positive | negative |
| Trial reagent | positive | a | b |
| | negative | c | d |

Calculation with the following formula:

Positive percent agreement = $a / (a+c) \times 100\%$

Negative percent agreement = $d / (b+d) \times 100\%$

Total percent agreement = $(a+d) / (a+b+c+d) \times 100\%$

5.2 Calculation method of 95% confidence interval

(1) The formula for calculating 95% confidence interval of total coincidence rate

$$[100\% (Q1 - Q2) / Q3, 100\% (Q1 + Q2) / Q3]$$

$$Q1 = 2(a+d) + 1.96^2$$

$$Q2 = 1.96 \sqrt{1.96^2 + 4(a+d)(b+c)/n}$$

$$Q3 = 2(n + 1.96^2)$$

(2) The formula for calculating 95% confidence interval of PPA

$$[100\% (Q1, ppa - Q2, ppa) / Q3, 100\% (Q1, ppa + Q2, ppa) / Q3, ppa]$$

$$Q1, ppa = 2a + 1.96^2$$

$$Q2, ppa = 1.96 \sqrt{1.96^2 + 4ac/(a+c)}$$

$$Q3, ppa = 2(a+c + 1.96^2)$$

(3) The formula for calculating 95% confidence interval of NPA

$$[100\% * (Q1, npa - Q2, npa) / Q3, npa, 100\% * (Q1, npa + Q2, npa) / Q3, npa]$$

$$Q1, npa = 2d + 1.96^2$$

$$Q2, npa = 1.96 \sqrt{1.96^2 + 4bd/(b+d)}$$

$$Q3, npa = 2(b+d + 1.96^2)$$

6. Test results

6.1 Clinical agreement study result

Across all study sites, a total of 257 healthy people samples and 152 volunteer's samples in serum collected at 28 days after the second dose of vaccine were tested with SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) and one comparator reagent (cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit). Overall study results are shown in Table 4 below.

Table 4 The test result of clinical agreement study

| Sample ID | GenScript | Getein | Sample ID | GenScript | Getein |
|-----------|-----------|---------|-----------|-----------|---------|
| S502151 | 58% | 67.4% | S472070 | 1% | < 10.0% |
| S473088 | 3% | < 10.0% | S473364 | 14% | 25.8% |
| S501627 | 47% | 58.2% | S475313 | 0% | < 10.0% |
| S564542 | 0% | < 10.0% | S567357 | 52% | 58.7% |
| S506297 | 76% | 84.8% | S502824 | 36% | 37.6% |
| S479002 | 5% | 12.6% | S470007 | 1% | < 10.0% |
| S547035 | 73% | 74.6% | S474360 | 4% | < 10.0% |
| S477572 | 6% | < 10.0% | S562458 | 82% | 80.8% |
| S477547 | 1% | < 10.0% | S472829 | 9% | < 10.0% |
| S476683 | 7% | 25.6% | S477760 | 17% | < 10.0% |
| S471044 | 13% | 22.3% | S474936 | 18% | < 10.0% |
| S564136 | 56% | 52.2% | S511162 | 37% | 32.9% |
| S563778 | 80% | 81.8% | S477759 | 6% | < 10.0% |
| S471613 | 18% | < 10.0% | S470870 | 0% | < 10.0% |
| S478008 | 7% | < 10.0% | S559789 | 55% | 69.3% |
| S472016 | 5% | 12.7% | S472104 | 4% | < 10.0% |
| S476990 | 0% | < 10.0% | S557538 | 78% | 76.1% |
| S545710 | 33% | 38.0% | S473087 | 11% | < 10.0% |
| S471969 | 4% | < 10.0% | S478071 | 11% | < 10.0% |
| S526734 | 15% | 12.6% | S560772 | 60% | 71.5% |
| S557509 | 61% | 58.7% | S472163 | 5% | 34.6% |
| S505999 | 36% | 36.9% | S475595 | 0% | < 10.0% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S553136 | 39% | 47.7% | S473218 | 0% | 15.6% |
| S569672 | 48% | 57.9% | S470629 | 10% | 11.6% |
| S473643 | 0% | < 10.0% | S541180 | 6% | < 10.0% |
| S470403 | 8% | < 10.0% | S473500 | 16% | < 10.0% |
| S479875 | 1% | < 10.0% | S557341 | 61% | 54.3% |
| S478650 | 0% | < 10.0% | S474568 | 13% | < 10.0% |
| S475873 | 14% | < 10.0% | S475451 | 7% | < 10.0% |
| S543837 | 13% | 21.7% | S511206 | 55% | 57.2% |
| S478014 | 0% | < 10.0% | S475938 | 15% | < 10.0% |
| S471417 | 1% | 20.1% | S531101 | 78% | 69.7% |
| S471034 | 9% | < 10.0% | S477793 | 11% | < 10.0% |
| S475023 | 18% | < 10.0% | S523997 | 37% | 37.5% |
| S474195 | 14% | 15.5% | S470387 | 0% | 12.6% |
| S527532 | 69% | 60.4% | S524045 | 89% | 81.6% |
| S470593 | 12% | < 10.0% | S479141 | 0% | < 10.0% |
| S471619 | 12% | < 10.0% | S478618 | 9% | 14.6% |
| S473104 | 0% | < 10.0% | S567996 | 4% | < 10.0% |
| S506292 | 87% | 90.9% | S479223 | 0% | < 10.0% |
| S477558 | 0% | < 10.0% | S475617 | 7% | < 10.0% |
| S472570 | 7% | < 10.0% | S477336 | 8% | < 10.0% |
| S478990 | 8% | 24.4% | S536878 | 41% | 41.6% |
| S478549 | 0% | < 10.0% | S470583 | 12% | 22.3% |
| S473752 | 17% | 11.2% | S507842 | 45% | 49.5% |
| S470368 | 10% | < 10.0% | S512594 | 40% | 49.8% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S476196 | 8% | 18.7% | S537643 | 45% | 25.9% |
| S539696 | 68% | 83.1% | S473240 | 11% | 14.3% |
| S557476 | 41% | 37.8% | S473374 | 5% | < 10.0% |
| S473817 | 0% | 21.8% | S478448 | 7% | 19.7% |
| S475919 | 18% | 19.2% | S518685 | 64% | 59.0% |
| S513164 | 63% | 49.3% | S533450 | 71% | 59.2% |
| S539063 | 57% | 68.4% | S477744 | 16% | 19.4% |
| S504628 | 45% | 41.1% | S471304 | 4% | < 10.0% |
| S479394 | 8% | 16.2% | S557798 | 61% | 61.3% |
| S550233 | 47% | 41.9% | S560602 | 80% | 86.2% |
| S472855 | 16% | 12.3% | S500909 | 55% | 44.2% |
| S477625 | 4% | < 10.0% | S470143 | 15% | < 10.0% |
| S470556 | 11% | < 10.0% | S500522 | 52% | 63.6% |
| S530823 | 83% | 77.8% | S564978 | 57% | 49.5% |
| S557460 | 90% | 93.5% | S520068 | 47% | 57.9% |
| S530763 | 88% | 86.7% | S474983 | 12% | 16.7% |
| S473024 | 5% | 14.0% | S472752 | 2% | < 10.0% |
| S475657 | 7% | < 10.0% | S520230 | 68% | 86.3% |
| S565421 | 74% | 85.1% | S510532 | 0% | < 10.0% |
| S473021 | 16% | < 10.0% | S476254 | 1% | < 10.0% |
| S472433 | 14% | < 10.0% | S477935 | 18% | < 10.0% |
| S475500 | 6% | < 10.0% | S479962 | 7% | 16.3% |
| S476554 | 0% | < 10.0% | S506056 | 35% | 20.1% |
| S476119 | 4% | < 10.0% | S534916 | 59% | 64.4% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S551231 | 69% | 84.2% | S470025 | 0% | 16.3% |
| S470900 | 1% | < 10.0% | S470924 | 5% | < 10.0% |
| S471025 | 14% | < 10.0% | S563651 | 71% | 57.7% |
| S516760 | 7% | < 10.0% | S564448 | 78% | 96.7% |
| S473722 | 0% | < 10.0% | S557285 | 65% | 73.2% |
| S550045 | 52% | 47.5% | S503466 | 57% | 68.2% |
| S471872 | 7% | < 10.0% | S472531 | 4% | < 10.0% |
| S475104 | 0% | < 10.0% | S476201 | 11% | < 10.0% |
| S473573 | 12% | < 10.0% | S477376 | 4% | 10.8% |
| S474914 | 12% | < 10.0% | S478734 | 19% | < 10.0% |
| S511845 | 7% | < 10.0% | S475990 | 9% | 26.3% |
| S470785 | 14% | < 10.0% | S547249 | 60% | 67.3% |
| S477940 | 12% | < 10.0% | S531674 | 38% | 46.6% |
| S555694 | 73% | 65.3% | S470355 | 3% | < 10.0% |
| S474178 | 0% | < 10.0% | S541848 | 58% | 61.9% |
| S523242 | 58% | 69.8% | S477221 | 6% | < 10.0% |
| S470008 | 10% | 10.4% | S475161 | 1% | < 10.0% |
| S514023 | 93% | 82.6% | S518563 | 52% | 46.1% |
| S532254 | 47% | 47.8% | S478046 | 15% | < 10.0% |
| S476940 | 15% | 15.2% | S552794 | 60% | 73.9% |
| S524989 | 57% | 49.5% | S470609 | 14% | < 10.0% |
| S476956 | 15% | < 10.0% | S476386 | 1% | < 10.0% |
| S540914 | 88% | 89.6% | S500910 | 47% | 45.7% |
| S477613 | 3% | < 10.0% | S514071 | 68% | 71.5% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S473156 | 6% | 12.1% | S471552 | 7% | < 10.0% |
| S533744 | 47% | 51.8% | S565364 | 41% | 40.6% |
| S477739 | 18% | < 10.0% | S478397 | 9% | < 10.0% |
| S477205 | 8% | < 10.0% | S473338 | 4% | 19.4% |
| S477837 | 12% | < 10.0% | S473257 | 2% | 12.9% |
| S479186 | 3% | 24.4% | S475542 | 4% | < 10.0% |
| S511159 | 92% | 95.9% | S476986 | 1% | < 10.0% |
| S479521 | 17% | < 10.0% | S478389 | 15% | 26.0% |
| S554527 | 88% | 83.0% | S478393 | 0% | < 10.0% |
| S557566 | 52% | 59.2% | S470075 | 18% | < 10.0% |
| S477591 | 3% | 22.3% | S471124 | 1% | 15.2% |
| S475000 | 11% | < 10.0% | S470807 | 1% | < 10.0% |
| S476163 | 14% | 11.6% | S540804 | 15% | 14.3% |
| S555978 | 75% | 93.0% | S567922 | 80% | 86.9% |
| S477568 | 0% | < 10.0% | S533004 | 66% | 59.8% |
| S479693 | 11% | < 10.0% | S565909 | 58% | 59.5% |
| S474535 | 9% | < 10.0% | S470056 | 7% | 15.3% |
| S551319 | 95% | 75.7% | S474966 | 13% | < 10.0% |
| S472826 | 1% | 16.2% | S479133 | 6% | 21.3% |
| S479254 | 3% | < 10.0% | S552009 | 34% | 36.1% |
| S479566 | 5% | < 10.0% | S475388 | 2% | < 10.0% |
| S527507 | 49% | 46.9% | S479106 | 12% | < 10.0% |
| S472941 | 10% | < 10.0% | S476955 | 0% | 14.9% |
| S516012 | 14% | < 10.0% | S476234 | 0% | < 10.0% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S472647 | 14% | < 10.0% | S566038 | 49% | 47.0% |
| S470313 | 4% | < 10.0% | S476871 | 16% | < 10.0% |
| S521281 | 77% | 84.4% | S476304 | 5% | < 10.0% |
| S479392 | 0% | 14.6% | S472337 | 18% | 17.0% |
| S479668 | 4% | < 10.0% | S475643 | 9% | < 10.0% |
| S473027 | 17% | < 10.0% | S477686 | 14% | 26.6% |
| S532124 | 68% | 63.2% | S471681 | 0% | < 10.0% |
| S478093 | 1% | 13.8% | S505230 | 21% | < 10.0% |
| S552736 | 93% | 79.5% | S568011 | 66% | 73.0% |
| S479357 | 0% | 15.3% | S472751 | 2% | < 10.0% |
| S476584 | 16% | < 10.0% | S556215 | 49% | 51.4% |
| S501436 | 11% | < 10.0% | S522183 | 45% | 39.9% |
| S479487 | 3% | < 10.0% | S471383 | 0% | < 10.0% |
| S477258 | 0% | 13.9% | S474765 | 12% | < 10.0% |
| S473558 | 12% | 11.30% | S473639 | 0% | < 10.0% |
| S473112 | 0% | < 10.0% | S472510 | 8% | < 10.0% |
| S477997 | 1% | 10.4% | S567144 | 91% | 83.4% |
| S548614 | 46% | 52.9% | S577513 | 0% | < 10.0% |
| S470698 | 17% | < 10.0% | S535409 | 0% | 18.8% |
| S478444 | 11% | 12.6% | S477800 | 0% | < 10.0% |
| S510608 | 66% | 73.6% | S478073 | 2% | 16.4% |
| S475810 | 5% | < 10.0% | S471601 | 10% | 22.3% |
| S470739 | 2% | 15.3% | S475042 | 19% | < 10.0% |
| S535624 | 65% | 61.5% | S473196 | 0% | < 10.0% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S477513 | 17% | < 10.0% | S552563 | 53% | 50.6% |
| S555403 | 69% | 60.7% | S470621 | 6% | 12.3% |
| S476855 | 7% | < 10.0% | S539775 | 57% | 55.0% |
| S477819 | 16% | < 10.0% | S530458 | 76% | 67.8% |
| S475946 | 11% | < 10.0% | S539594 | 35% | 37.2% |
| S477701 | 2% | < 10.0% | S477991 | 0% | < 10.0% |
| S479358 | 4% | < 10.0% | S475556 | 8% | < 10.0% |
| S508357 | 60% | 63.3% | S477938 | 19% | 20.0% |
| S474417 | 13% | 23.0% | S479543 | 17% | < 10.0% |
| S477504 | 1% | 14.4% | S562071 | 23% | 11.2% |
| S470163 | 7% | 19.7% | S510089 | 34% | 32.9% |
| S535245 | 45% | 39.2% | S514157 | 59% | 50.2% |
| S475735 | 16% | < 10.0% | S478189 | 10% | < 10.0% |
| S554291 | 70% | 76.4% | S470750 | 18% | < 10.0% |
| S472399 | 11% | 20.1% | S476354 | 14% | < 10.0% |
| S547257 | 84% | 78.4% | S474886 | 3% | < 10.0% |
| S475722 | 9% | < 10.0% | S471641 | 0% | < 10.0% |
| S478069 | 8% | < 10.0% | S479124 | 14% | 20.1% |
| S477199 | 9% | < 10.0% | S478298 | 6% | < 10.0% |
| S542380 | 37% | 39.5% | S472566 | 12% | 25.0% |
| S473915 | 9% | 19.2% | S474724 | 10% | < 10.0% |
| S477965 | 18% | < 10.0% | S538188 | 65% | 60.0% |
| S474991 | 4% | < 10.0% | S474205 | 7% | 20.7% |
| S553199 | 11% | 16.2% | S545502 | 47% | 43.0% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S506840 | 95% | 90.0% | S519000 | 93% | 78.1% |
| S507353 | 65% | 54.3% | S534609 | 44% | 38.7% |
| S473414 | 10% | < 10.0% | S509082 | 56% | 59.1% |
| S472151 | 16% | < 10.0% | S474369 | 14% | 15.4% |
| S476116 | 15% | < 10.0% | S520093 | 5% | < 10.0% |
| S475687 | 17% | < 10.0% | S479578 | 6% | 21.5% |
| S476493 | 18% | < 10.0% | S470290 | 6% | < 10.0% |
| S477450 | 13% | 15.0% | S514891 | 58% | 61.9% |
| S472230 | 18% | < 10.0% | S532216 | 87% | 87.2% |
| S555117 | 61% | 57.9% | S527460 | 5% | < 10.0% |
| S472559 | 16% | < 10.0% | S473474 | 15% | < 10.0% |
| S473363 | 2% | < 10.0% | S477639 | 6% | 15.0% |
| S471033 | 7% | 12.6% | S511158 | 80% | 77.4% |
| S479780 | 6% | < 10.0% | S472397 | 0% | 25.6% |
| S520221 | 54% | 62.0% | S500710 | 11% | < 10.0% |
| S470598 | 16% | 10.80% | S511877 | 13% | < 10.0% |
| S470125 | 5% | 13.1% | S471454 | 17% | < 10.0% |
| S501290 | 89% | 83.5% | S509345 | 3% | 11.9% |
| S472014 | 7% | < 10.0% | S471263 | 6% | 11.1% |
| S549251 | 38% | 41.5% | S475208 | 4% | < 10.0% |
| S470273 | 6% | < 10.0% | S472157 | 0% | < 10.0% |
| S474737 | 13% | < 10.0% | S471024 | 0% | 21.2% |
| S475911 | 1% | < 10.0% | S521406 | 57% | 58.2% |
| S530192 | 67% | 63.6% | S550714 | 43% | 37.6% |

| | | | | | |
|---------|-----|---------|---------|-----|---------|
| S476410 | 11% | < 10.0% | S476663 | 11% | 20.2% |
| S516304 | 60% | 54.8% | S513867 | 85% | 83.8% |
| S479257 | 12% | < 10.0% | S476291 | 5% | 26.1% |
| S476318 | 15% | < 10.0% | S479301 | 5% | < 10.0% |
| S555810 | 66% | 56.6% | S479113 | 0% | < 10.0% |
| S568154 | 82% | 80.6% | S478292 | 13% | < 10.0% |
| S471125 | 16% | < 10.0% | S476713 | 16% | < 10.0% |
| S470447 | 13% | 18.0% | S476392 | 1% | < 10.0% |
| S567823 | 82% | 90.4% | S472240 | 0% | < 10.0% |
| S472317 | 14% | < 10.0% | S476907 | 0% | 12.6% |
| S529420 | 34% | 35.5% | S471040 | 12% | < 10.0% |
| S475576 | 9% | 12.6% | S542741 | 53% | 51.5% |
| S472975 | 0% | 17.4% | S542825 | 46% | 38.9% |
| S558088 | 81% | 77.2% | S475441 | 13% | < 10.0% |
| S557506 | 87% | 98.2% | | | |

Note: For GenScript Biotech Co., Ltd cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit, the test result is displayed numerically in terms of inhibition ratio. Test result is negative if inhibition rate is <30%, and positive if inhibition rate is ≥ 30%.

For Getein Biotech, Inc. SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay), Test result is negative if inhibition rate is <30%, and positive if inhibition rate is ≥ 30%.

7. Statistical analysis

7.1 Clinical agreement study results analysis

According to the statistical methods, the results of clinical agreement study were as follows.

Correlation analysis of SARS-CoV-2 neutralizing antibody

The test results were summarized as follows.

Table 5 Overall clinical study results for SARS-CoV-2 neutralizing antibody

| | |
|----------------------------------|--------------------------------|
| SARS-CoV-2 Neutralizing Antibody | Comparator reagent (GenScript) |
|----------------------------------|--------------------------------|

| | | positive | negative |
|---------------|----------|----------|----------|
| Trail reagent | positive | 130 | 0 |
| | negative | 1 | 278 |

Compared to GenScript neutralizing antibody (cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit):

Positive Percent Agreement = $130 / (130+1) \times 100\% = 99.2\%$ (95%CI: 95.80% - 99.87%)

Negative Percent Agreement = $278 / (278+0) \times 100\% = 100.0\%$ (95%CI: 98.64% - 100.00%)

Total percent agreement = $(130+278) / (130+0+1+278) \times 100\% = 99.8\%$ (95%CI: 98.63% - 99.96%)

8. Conclusion

In this study, the trail reagent is SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay), developed by Getein Biotech, Inc. The the comparator reagent is cPass™ SARS-CoV-2 Neutralization Antibody Detection Kit *for detecting SARS-CoV-2* neutralizing antibodies developed by GenScript Biotech Co., Ltd. A total of 409 individual patient samples, with 257 healthy people and 152 volunteers who have been immunized by injecting two doses of vaccine against SARS-CoV-2, were tested with the trial reagent to evaluate the clinical agreement with the comparator.

According to the results of clinical agreement study, 409 clinical patients were tested. For SARS-CoV-2 Neutralizing Antibody tested by trial reagents, the positive percent agreement was 99.2%, the 95% confidence interval of the positive percent agreement was [95.80%,99.87%], the negative percent agreement was 100.0%, the 95% confidence interval of the negative percent agreement was [98.64%,100.00%], the total percent agreement is 99.8%, the 95% confidence interval of the total percent agreement was [98.63%,99.96%].